
OVA & PARASITE

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
World Health Organization. Bench aids for the diagnosis of intestinal parasites. WHO: Geneva, 1994.
Garcia LS, et al. Practical Guidance for Clinical Microbiology Laboratories: Laboratory Diagnosis of Parasites from the Gastrointestinal Tract. Clin Microbiol 2018 Rev 31:10.1128/cmr.00025-17.
Branda JA, et al. A rational approach to the stool ova and parasite examination. Clin Infect Dis. 2006 Apr 1;42(7):972-8.

OVA & PARASITE EXAMINATION

Kanyapak Kooburat
Clinical Pathologist, Pathology dept.

A routine ova and parasite examination (O&P) usually includes laboratory procedures that are designed to detect organism in clinical specimens by using microscopic characteristics rather than culture, biochemical tests, and physical growth characteristics.

The test that is to be requested depends on which infectious agent is suspected. It is very important that the ordering physician be aware of what procedure are performed with the O&P and that this test may not cover all parasites. Additional stains for Microsporida and *Cyclospora*/*Cystoisospora*/*Cryptosporidium* may need to be ordered separately, and immunoassays for selected organisms may also not be included in the O&P.

Patient preparation

- Should be performed **before barium is administered** to a patient for radiology exams.
- Stool specimens **containing the opaque, chalky suspension are unacceptable**, and intestinal protozoa may be masked for 5-10 days after ingestion of barium.
- **Caster oil or mineral oil, bismuth, antibiotic including antimalarial medication, and non-absorbable antidiarrheal preparations, may interfere** examination, and collection should be postponed for 5-10 days after administration to allow clearance of these substances.

Best examination time

Fresh stool specimens are mandatory of motile protozoan trophozoites (amoebae, flagellates, ciliates)

- **Liquid stools** should be examined within 30 min of passage (not after arriving in the lab).
- **Semi-formed or soft stools** should be examined within 1 hr of passage, usually contain both cyst and trophozoites.
- **Formed stools** should be examined within 24 hr after passage and contain mainly cysts.


Things must be known

****It is recommended that multiple stool samples be examined prior to ruling out a parasitic infection.**

******* Many organisms, particularly intestinal protozoa, do not appear in the stool in consistent numbers. **Physicians should be aware that the probability of detecting clinically relevant parasites in a single specimen may be low as 50-60% but is >95% if three samples are examined by O&P.**



What to do when: Delayed handling specimens

- **Storing stool at 4–8°C** keeps parasites alive for different lengths of time depending on their stage. *E. histolytica* trophozoites survive 1–3 days, cysts about 4 weeks, and helminth eggs stored at 3–5°C can remain viable for up to 1 month.

Direct wet-mount prep

- Prepared by mixing a small amount of fresh, unpreserved stool with saline then suspension under coverslip for motile protozoan trophozoites.
- **Helminth eggs, protozoan cysts, and *Cyclospora/Cystoisospora/Cryptosporidium* oocysts may also be observed on wet film.**
- Motility can not be observed from preserved specimens

Concentration procedure

Fecal conc. is a routine part of a complete O&P. It allows detection of small numbers of organisms that may be missed on a direct wet mount.

Two type of concentration procedure:

1. **Sedimentation (FEA):** by centrifugation parasites and exams with or without iodine
2. **Flotation (zinc sulfate):** used high specific gravity of the solution floats the organisms

Both are **used to concentrate helminth eggs and larvae, protozoan cysts, *Cyclospora/Cystoisospora/Cryptosporidium* oocysts, and microsporidial spores.**


EXAMINATION METHODS

The most commonly performed test in the parasitology laboratory is the complete O&P.

- Direct examination: gross specimen and physical characteristics can provide information. Worms may be seen in the stool specimen and retrieved for identification.
- **Microscopic examination:** Direct wet-mount preparation, Concentration procedures, Routine permanent-stain smears
- Modified trichrome stain
- Iron-hematoxylin

Other special methods for *Cyclospora*, *Cystoisospora*, *Cryptosporidium* such as Special stains, Modified Acid-Fast staining, Concentration of *Cryptosporidium* and *Cyclospora*, Safranin staining or Fluorescent antibody staining (DFA assay) for *Cryptosporidium*, Modified trichrome stains for chromotrope 2R.

Modified trichrome stain

The stain is easy to perform and allows detection of protozoan trophozoites and cysts. The color stain allowing easier differentiation between organisms and artifacts

Iron-hematoxylin

- Two most common used are Spencer-Monroe and Tompkins-Miller procedures
- **Used for most of the intestinal protozoans**
- Specimen stain shades of grayish blue
- Can be used in preservative specimens

Take Notes


- Do not order an O&P for the diagnosis of Enterobiasis.
- Pin worm prep: use to collect eggs from the perianal skin

References

World Health Organization. Bench aids for the diagnosis of intestinal parasites. WHO: Geneva, 1994.
Garcia LS, et al. Practical Guidance for Clinical Microbiology Laboratories: Laboratory Diagnosis of Parasites from the Gastrointestinal Tract. Clin Microbiol 2018 Rev 31:10.1128/cmr.00025-17.
Branda JA, et al. A rational approach to the stool ova and parasite examination. Clin Infect Dis. 2006 Apr 1;42(7):972-8.

OVA & PARASITE 检查

Kanyapak Kooburat
Clinical Pathologist, Pathology dept.

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
References
World Health Organization. Bench aids for the diagnosis of intestinal parasites. WHO: Geneva, 1994.
Garcia LS, et al. Practical Guidance for Clinical Microbiology Laboratories: Laboratory Diagnosis of Parasites from the Gastrointestinal Tract. Clin Microbiol 2018 Rev 31:10.1128/cmr.00025-17.
Branda JA, et al. A rational approach to the stool ova and parasite examination. Clin Infect Dis. 2006 Apr 1;42(7):972-8.